Belluscura LLC 8-cell Battery

UN 38.3 Lithium Battery Test Summary

Cell, Battery or Product Model Number				
Item Number:	50188			
Item Name:	8-cell Battery			

Cell, battery, or product manufacturer's contact information							
Name:	Belluscura LLC						
Address:	5504 Democracy Drive, Suite 200						
City:	Plano	State:	TX	ZIP:	75024	Country:	USA
Telephone:	(888) 469-7567	E-ma	ailcustomerse	ervice@bellus	scura.com/ebsit	e: https://www.bel	luscura.com/

Test Labo	ratory				
Name:	SGS				
Address:	620 Old Peachtree Road NW, Suite 100				
City:	Suwanee	State: GA	ZIP:	30024	Country: USA
Telephone:	(770) 570-1800	E-mail:	Website: https://www.sgs.com/		

Cell	or Ba	ttery Description					
Cell or Battery: X-PLO2R 8-cell					Physical Description (dimensions, appearance):		
Cell c	Cell or Battery Type: lithium-ion				18.5 cm x 7.6 cm x 2.9 cm		
Watt-	Watt-hour rating or Lithium Content: 99.36 watt-hours			ours	black plastic case charge port, LED display		
Completed Cell or Battery Weight: 562 grams			62 grams				
Uniqu	Unique Test Report ID Number: SUW-4544691-2			91-2	Date of test report: May 31, 2019		
		List of Tests Completed:					
Yes	No		Pass	Fail	Additional Comments (or indicate compliance with other standards, e.g., Underwriters Laboratory):		
Х		Test T.1: Altitude simulation	Х				
Х		Test T.2: Thermal test	Х				
Х		Test T.3: Vibration	Х		Reference to assembled battery testing requirements, if applicable (i.e., 38.3.3 (f) and 38.3.3 (g)).		
Х		Test T.4: Shock	Х		UN DOT 38.3		
Х		Test T.5: External short circu	uit X		IEC 62133-2:2017		
	Х	Test T.6: Impact/Crush			Reference to the revised edition of the Manual of Tests and Criteria used and to amendments thereto, if any:		
Х		Test T.7: Overcharge	Х		TRANSPORT OF DANGEROUS GOODS		
	Х	Test T.8: Forced discharge			Manual of Tests and Criteria, 6th Revised Edition, Amendment 1		

Signature with name and title of signatory as an indication of the validity of information provided: *R. Paul Bray*, Vice President, Quality Date:

March 2, 2022